

Products to Maximize Human Resources

HEMASTAIN PRODUCTS

1. SedonaLab.com's Hemastain

Our Flagship Hemastain (Wright's Stain) is the best staining solution available for peripheral blood films. Every lot is tested to determine absorbance characteristics at various wavelengths of light. The stain is then qualified on reference Hematrak systems to assure accuracy on all manual and automated microscopy, including imaging systems. This continual testing allows for lot to lot conformity. Uniquely effective when used in conjunction with our Hemastain Fixative and Phosphate Buffer.



2. Stain Fixative

Every lot of our Hemastain fixative is tested with our Hemastain (Wright's Stain) and Phosphate Buffer to ensure quality results. We believe our Stain Fixative is your best solution for uniform staining results.

3. Buffer Powder

The Phosphate Buffer is pre-measured and certified to meet our exacting standards. Easy to use, pre-weighed in foil packets to make one gallon of buffer solution. Our buffer powder is manufactured to maintain a specific pH throughout its shelf life (7.1pH), assuring optimum staining quality every time. Differentiation is further enhanced by a final buffer rinse, a SedonaLab.com exclusive process which provides a unique degree of quality and clarity in the staining process.

AUTOPREP®



Features:

- User Friendly
- Easy Setup, Use and Clean-up
- Constant Dispense Rate
- Standard Wedge Configuration
- Replaceable Glass Blades
- Low Maintenance
- Stationary Or Portable Use
- Smear/Speed Control Knob
- Pause Control
- Uses Fresh or Anti-Coagulated Blood
- Durable Construction, Stainless Steel Case
- Requires no Power Source or Battery
- Extra Glass Blades Included
- One Year Warranty

MINIPREP PLUS®



Step 1: Place slides




Step 2: Drop blood



Step 3: Press down



Step 4: Release



Step 5: Clean glass

Visit our online store at www.SedonaLab.com

2370 West Hwy. 89A, Suite 11-507, Sedona, AZ 86336

Phone and Fax: 1-800-851-9844

Email: info@SedonaLab.com